




# Fine Particle Filter Installation Instructions

en

Instructions For Use

<b>1</b> <b>➔ Remove the filter from the package.</b>	
<b>2</b> <b>➔ Place the filter into the housing with the foam (dark) side facing toward the ventilator.</b> <b>➔ Ensure that the filter is flush with the inside surface of the housing.</b>	
<b>3</b> <b>➔ Installation complete.</b>	

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.



Single use


Caution, consult  
accompanying  
documents

Part No. 10072519 Rev B 2014-07

 GVS Filter Technology UK Ltd  
 NFC House, Vickers Industrial Estate  
 Mellishaw Lane, Morecambe  
 Lancashire, LA3 3EN, United Kingdom  
 Distributed by Covidien Ilc 15 Hampshire Street  
 Mansfield, MA 02048 USA.

[www.covidien.com](http://www.covidien.com)